



Glowing products for science

Call us : [800-304-5357](tel:800-304-5357)

## CF® Dye SE/TFP Ester

CF® Dye SE/TFP esters are amine-reactive fluorescent dyes. They are commonly used to label primary amines of antibodies and other proteins, amine-modified oligonucleotides and other amine-containing molecules.



### Product attributes

<b>Chemical reactivity (reacts with)</b>	Amines
<b>Functional group</b>	Tetrafluorophenyl Ester (TFP Ester), Succinimidyl Ester (NHS Ester)
<b>Storage Conditions</b>	Store at -10 to -35 °C, Protect from light

## Product Description

CF® Dye SE/TFP Esters are amine-reactive forms of Biotium's bright and photostable CF® Dyes. The succinimidyl ester (SE) group, also known as N-hydroxysuccinimide (NHS) esters, on the dye molecule reacts with primary amines such as lysine residues of antibodies or proteins, amine-modified oligonucleotides and other amine-containing molecules to form a stable amide linkage. CF® Dye TFP (tetrafluorophenyl) esters are more stable alternatives to succinimidyl ester. CF® Dye TFP esters are available for CF®850 and CF®870, Biotium's industry-leading near-infrared dyes with emission above 850 nm.

- Amine-reactive fluorescent dyes with high labeling efficiency
- Label antibodies, proteins, amine-modified oligonucleotides and other amine-containing molecules.
- Offered with Biotium's bright and photostable CF® Dyes which provide superior signal-to-noise
- Available in over 40 CF® Dye colors from blue to near-IR wavelengths
- *Near-IR CF®850 and CF®870 TFP esters offer emission above 850 nm*

Due to Biotium's unique manufacturing and packaging process, our amine-reactive CF® Dyes are provided with the highest purity and more stable ester dye forms than other vendors. As a result, our CF® Dye SE/TFP products retain more reactivity and provide a higher labeling efficiency for users.

### Superior CF® Dyes

CF® Dyes are Biotium's line of next-generation fluorescent dyes that have improved brightness, photostability, and water solubility compared to other commercially available fluorescent dyes. Learn more about [CF® Dyes](#). For more information download the [CF® Dye Brochure](#).

Also see our [CF® Dye SE Protein Labeling Kits](#) and [CF® Dye Mix-n-Stain™ Antibody Labeling Kits](#).

## CF® Dye SE/TFP Ester

Dye	Ex/Em	Size	Catalog No.	Dye Features
<a href="#">CF@350 SE</a>	347/448 nm	1 umol	<a href="#">92109</a>	<a href="#">CF@350 Features</a>
<a href="#">CF@405S SE</a>	404/431 nm	1 umol	<a href="#">92110</a>	<a href="#">CF@405S Features</a>
<a href="#">CF@405M SE</a>	408/452 nm	1 umol	<a href="#">92111</a>	<a href="#">CF@405M Features</a>
<a href="#">CF@405L SE</a>	395/545 nm	1 umol	<a href="#">92112</a>	<a href="#">CF@405L Features</a>
<a href="#">CF@410 SE</a>	404/455 nm	5 mg	<a href="#">97502</a>	
<a href="#">CF@430 SE</a>	426/498 nm	1 umol	<a href="#">92117</a>	<a href="#">CF@430 Features</a>
<a href="#">CF@440 SE</a>	440/515 nm	1 umol	<a href="#">92123</a>	<a href="#">CF@440 Features</a>
<a href="#">CF@450 SE</a>	450/538 nm	1 umol	<a href="#">96011</a>	<a href="#">CF@450 Features</a>
<a href="#">CF@488A SE</a>	490/515 nm	1 umol	<a href="#">92120</a>	<a href="#">CF@488A Features</a>
<a href="#">CF@503R SE</a>	503/532 nm	1 umol	<a href="#">96078</a>	<a href="#">CF@503R Features</a>
<a href="#">CF@505 SE</a>	505/519 nm	1 umol	<a href="#">97500</a>	Structurally equivalent to ATTO 488
<a href="#">CF@510 SE</a>	513/537 nm	1 umol	<a href="#">96146</a>	
<a href="#">CF@514 SE</a>	516/548 nm	1 umol	<a href="#">92103</a>	<a href="#">CF@514 Features</a>
<a href="#">CF@532 SE</a>	527/558 nm	1 umol	<a href="#">92104</a>	<a href="#">CF@532 Features</a>
<a href="#">CF@543 SE</a>	541/560 nm	1 umol	<a href="#">92105</a>	<a href="#">CF@543 Features</a>
<a href="#">CF@550R SE</a>	551/577 nm	1 umol	<a href="#">96073</a>	<a href="#">CF@550R Features</a>
<a href="#">CF@555 SE</a>	555/565 nm	1 umol	<a href="#">92130</a>	<a href="#">CF@555 Features</a>
<a href="#">CF@568 SE</a>	562/583 nm	1 umol	<a href="#">92131</a>	<a href="#">CF@568 Features</a>
<a href="#">CF@570 SE</a>	568/591 nm	1 umol	<a href="#">96014</a>	<a href="#">CF@570 Features</a>
<a href="#">CF@583 SE</a>	583/606 nm	1 umol	<a href="#">96016</a>	<a href="#">CF@583 Features</a>
<a href="#">CF@583R SE</a>	586/609 nm	1 umol	<a href="#">96084</a>	<a href="#">CF@583R Features</a>
<a href="#">CF@594 SE</a>	593/614 nm	1 umol	<a href="#">92132</a>	<a href="#">CF@594 Features</a>
<a href="#">CF@597R SE</a>	597/619 nm	1 umol	<a href="#">96092</a>	<a href="#">CF@597R Features</a>
<a href="#">CF@620R SE</a>	617/639 nm	1 umol	<a href="#">92106</a>	<a href="#">CF@620R Features</a>
<a href="#">CF@633 SE</a>	630/650 nm	1 umol	<a href="#">92133</a>	<a href="#">CF@633 Features</a>
<a href="#">CF@640R SE</a>	642/662 nm	1 umol	<a href="#">92108</a>	<a href="#">CF@640R Features</a>
<a href="#">CF@647 SE</a>	650/665 nm	1 umol	<a href="#">92135</a>	<a href="#">CF@647 Features</a>
<a href="#">CF@660R SE</a>	663/682 nm	1 umol	<a href="#">92134</a>	<a href="#">CF@660R Features</a>
<a href="#">CF@660C SE</a>	667/685 nm	1 umol	<a href="#">92137</a>	<a href="#">CF@660C Features</a>
<a href="#">CF@680 SE</a>	681/698 nm	1 umol	<a href="#">92139</a>	<a href="#">CF@680 Features</a>
<a href="#">CF@680R SE</a>	680/701 nm	1 umol	<a href="#">92107</a>	<a href="#">CF@680R Features</a>
<a href="#">CF@700 SE</a>	695/720 nm	1 umol	<a href="#">96067</a>	<a href="#">CF@700 Features</a>
<a href="#">CF@710 SE</a>	712/736 nm	1 umol	<a href="#">96125</a>	<a href="#">CF@710 Features</a>
<a href="#">CF@725 SE</a>	729/750 nm	1 umol	<a href="#">96126</a>	<a href="#">CF@725 Features</a>
<a href="#">CF@740 SE</a>	742/767 nm	1 umol	<a href="#">96109</a>	<a href="#">CF@740 Features</a>
<a href="#">CF@750 SE</a>	755/779 nm	1 umol	<a href="#">92142</a>	<a href="#">CF@750 Features</a>
<a href="#">CF@770 SE</a>	770/797 nm	1 umol	<a href="#">92150</a>	<a href="#">CF@770 Features</a>
<a href="#">CF@790 SE</a>	784/806 nm	0.25 umol	<a href="#">92155</a>	<a href="#">CF@790 Features</a>
<a href="#">CF@800 SE</a>	797/816 nm	0.25 umol	<a href="#">92127</a>	<a href="#">CF@800 Features</a>
<a href="#">CF@820 SE</a>	822/835 nm	0.25 umol	<a href="#">96068</a>	<a href="#">CF@820 Features</a>
<a href="#">CF@850 TFP</a>	852/870 nm	0.25 umol	<a href="#">96094</a>	<a href="#">CF@850 Features</a>
<a href="#">CF@870 TFP</a>	876/896 nm	0.25 umol	<a href="#">96095</a>	<a href="#">CF@870 Features</a>

## CF® Dye Technical Data

Dye	Abs/Em (nm)	MW (free acid form)	Sephadex® media <sup>1</sup>	Amax (nm)	Cf A260/Amax	Cf A280/Amax	ε <sup>2</sup>	Optimal DOL (IgG)
<a href="#">CF@350</a>	347/448	~496	G-25	347	0.13	0.14	18,000	4-6
<a href="#">CF@405S</a>	404/431	~1169	G-25	404	0.19	0.7	33,000	5-10
<a href="#">CF@405M</a>	408/452	~503	G-25	408	0.24	0.13	41,000	4-6
<a href="#">CF@405L</a>	395/545	~1573	G-25	395	N/A	0.5	24,000	8-12
<a href="#">CF@410</a>	404/455	~242	G-25	416	0.15	0.2	46,000	5-7
<a href="#">CF@430</a>	426/498	~429	G-25	426	0.21	0.044	40,000	5-8
<a href="#">CF@440</a>	440/515	~479	G-25	440	0.303	0.139	40,000	5-8
<a href="#">CF@450</a>	405/460	~689	G-25	450	0.205	0.2	40,000	5-8
<a href="#">CF@488A</a>	490/515	~914	G-25	490	0.16	0.1	70,000	7-9
<a href="#">CF@503R</a>	503/532	~1100	G-25	503	0.21	0.09	90,000	4-10
<a href="#">CF@505</a>	505/519	~587	G-25	505	0.22	0.09	90,000	4-8
<a href="#">CF@510</a>	513/537	~562	G-25	513	0.25	0.14	90,000	3-5
<a href="#">CF@514</a>	516/548	~1216	G-25	516	0.14	0.073	105,000	5-8
<a href="#">CF@532</a>	527/558	~685	G-25	527	0.11	0.06	96,000	4-7
<a href="#">CF@543</a>	541/560	~887	G-25	541	0.305	0.095	100,000	4-7
<a href="#">CF@550R</a>	551/577	~686	G-25	551	0.12	0.08	100,000	5-6
<a href="#">CF@555</a>	555/565	~959	G-25	555	0.026	0.08	150,000	4-5, 3-6*

Dye	Abs/Em (nm)	MW (free acid form)	Sephadex® media <sup>1</sup>	Amax (nm)	Cf A260/Amax	Cf A280/Amax	ε <sup>2</sup>	Optimal DOL (IgG)
<a href="#">CF@568</a>	562/583	~714	G-25	562	0.17	0.08	100,000	5-6
<a href="#">CF@570</a>	568/591	~2998	G-25	568	0.0998	0.1	150,000	5-6
<a href="#">CF@583</a>	583/606	~3127	G-75	583	0.139	0.223	150,000	5-6
<a href="#">CF@583R</a>	585/609	~773	G-25	585	0.19	0.08	100,000	5-6
<a href="#">CF@594</a>	593/614	~729	G-25	593	0.24	0.08	115,000	4-7
<a href="#">CF@597R</a>	597/619	~800	G-25	597	0.25	0.08	100,000	5-6
<a href="#">CF@620R</a>	617/639	~738	G-25	617	0.28	0.45	115,000	5-6
<a href="#">CF@633</a>	630/650	~821	G-25	630	0.25	0.48	100,000	4-7
<a href="#">CF@640R</a>	642/662	~832	G-50	642	0.23	0.44	105,000	4-7
<a href="#">CF@647</a>	650/665	~985	G-25	650	0.01	0.03	240,000	4-5, 3-6*
<a href="#">CF@660C</a>	667/685	~3024	G-75	667	0.05	0.08	200,000	3-6, 2-3*
<a href="#">CF@660R</a>	663/682	~888	G-25	663	0.2	0.51	100,000	4-7
<a href="#">CF@680</a>	681/698	~3153	G-75	681	0.06	0.09	210,000	3-5, 2-3*
<a href="#">CF@680R</a>	680/701	~912	G-25	680	0.155	0.32	140,000	5-6
<a href="#">CF@700</a>	696/721	~2474	G-75	696	0.055	0.06	240,000	3-6
<a href="#">CF@710</a>	712/736	~860	G-25	712	0.11	0.07	115,000	5-6
<a href="#">CF@725</a>	729/750	~890	G-25	729	0.11	0.07	120,000	5-6
<a href="#">CF@740</a>	742/767	~900	G-25	742	0.132	0.08	105,000	5-6
<a href="#">CF@750</a>	755/777	~2921	G-75	755	0.01	0.03	250,000	3-5, 2-3*
<a href="#">CF@770</a>	770/797	~3091	G-75	770	0.041	0.06	220,000	3-5, 2-3*
<a href="#">CF@790</a>	784/806	~3179	G-75	784	0.104	0.07	210,000	3-5
<a href="#">CF@800</a>	797/816	~3334	G-75	797	0.09	0.08	210,000	3-5
<a href="#">CF@820</a>	822/835	~2711	G-75	822	0.0459	0.07	253,000	3-6
<a href="#">CF@850</a>	852/870	~2787	G-75	852	N/A	0.06	240,000	3-6
<a href="#">CF@870</a>	876/896	~2773	G-75	877	N/A	0.06	240,000	3-6

1. Sephadex® media recommendations are for antibody purification, not nucleic acid.

2. Extinction Coefficient (ε).

\*Suitable, but suboptimal DOL.

Sephadex is a registered trademark of Cytiva.

## References

1. Biomater. Sci.(2019) 7, 1898-1904. [DOI:10.1039/C8BM01378K](https://doi.org/10.1039/C8BM01378K)
2. Biomaterials (2015) 67, 382-392. [DOI:10.1016/j.biomaterials.2015.07.046](https://doi.org/10.1016/j.biomaterials.2015.07.046)
3. Microbiological Research (2018) 217, 69-80. [DOI:10.1016/j.micres.2018.08.017](https://doi.org/10.1016/j.micres.2018.08.017)

Download a list of [CF@ Dye references](#).

This datasheet was generated on May 8, 2026 at 03:06:07 PM. Visit product page to check for updated information before use.

Product link: <https://biotium-woo.supremeclients.com/product/cf-dye-succinimidyl-ester/>