

## Ethidium Monoazide Bromide (EMA)

A photo-activated fluorescent nucleic acid binding dye used to selectively and covalently label dead cells in the presence of live cells for viability PCR, fluorescent microscopy, and flow cytometry.



### Product Description

When not activated or bound to DNA, Ethidium monoazide bromide (EMA) is a non-fluorescent nucleic acid stain with a photoaffinity label. Once photoactivated, the dye binds covalently to nucleic acids and can emit fluorescence. The dye is most often used to label dead cells in a mixed population of live and dead cells.

Because EMA is relatively impermeant to live cells, it selectively binds DNA in dead cells. Photolysis renders the dead cell DNA covalently labeled with the dye. One can then wash and fix the cell preparation for analysis by microscopy, fluorescence plate reader or flow cytometry. The dye has been used to footprint drug binding sites on DNA, to modify plasmid DNA, and for labeling yeast cells for leukocyte phagocytosis assays. The major advantage of this method is that researchers can avoid extensive manipulation of live pathogenic organisms. EMA has been used to differentiate between viable and dead bacteria by 5'-nuclease PCR.

### Viability PCR

Switching from agar plates to v-PCR is a great way to speed up bacteria viability testing. There are some things that you will need to know and take into consideration before you switch to the v-PCR method. For viability PCR applications, [Ethidium monoazide \(EMA\)](#) is the most cell-permeable, and may result in false negatives. [Propidium monoazide \(PMA\)](#), developed at Biotium, is the most widely-used dye, so there are many protocols and publications available for reference. [PMAxx™](#) is designed to be superior to PMA and EMA as well as compatible with existing PMA protocols. [PMAxx™](#) has better differentiation of live and dead bacteria and therefore it is more selective for dead cell labeling for v-PCR quantitation of bacteria. See the [PMA & PMAxx™ Reference List](#) and [PMA & PMAxx™ Validated Bacterial Strains](#). For more information about v-PCR, [see our blog post](#).

### Viability PCR Kits

To more easily get started with viability PCR, learn about our customizable [Viability PCR Starter Kits](#). You choose which viability PCR dye you prefer, and whether or not you need Enhancer for gram-negative bacteria. All of the starter kits come with an EvaGreen® qPCR master mix and a detail protocol to get you started.

Biotium also provides strain-specific [Real-Time PCR Bacterial Viability Kits](#) with validated primers for selected pathogens: M. tuberculosis, S. aureus, MRSA, Salmonella, E. coli, E. coli O157:H7, Legionella and Listeria. These kits provide everything that you need for the selective detection of your favorite species of live bacteria by real-time v-PCR.

### LED Photolysis Devices

Biotium sells two devices for the photolysis step of viability PCR. The [PMA-Lite™ LED photolysis device](#) is designed for PMAxx™-, PMA-, and EMA-treated samples in microcentrifuge tubes. The [Glo-Plate™ Blue](#) is a flat LED box ideal for samples in microplates, filters, or other formats.

Materials from Biotium are sold for research use only.

### References

1. Plant Dis, 104, 1105, (2020), [DOI: 10.1094/PDIS-10-19-2248-RE](#)
2. Environ. Sci.: Water Res. Technol., 6, 87, (2020), [DOI: 10.1039/C9EW00896A](#)
3. J Leukoc Biol., Pre-Print, (2020), [DOI: 10.1002/JLB.5MA1219-296R](#)
4. Cell, 178, 1231, (2019), [DOI: 10.1016/j.cell.2019.07.033](#)
5. Curr Microbiol, 76, 1425, (2019), [DOI: 10.1007/s00284-019-01772-y](#)
6. Environ. Sci.: Water Res. Technol., 5, 91, (2019), [DOI: 10.1039/C8EW00613J](#)
7. AIMS Microbiol., 5, 102, (2019), [DOI: 10.3934/microbiol.2019.1.102](#)
8. Water, 11, 1196, (2019), [DOI: 10.3390/w11061196](#)
9. J Microbiol Methods, 148, 33, (2018), [DOI: 10.1016/j.mimet.2018.03.018](#)
10. Water, 10, 7, (2018), [DOI: 10.3390/w10070841](#)
11. Food control, 85, 23, (2018), [DOI: 10.1016/j.foodcont.2017.09.013](#)
12. Vaccine, 35, 4330, (2017), [DOI: 10.1016/j.vaccine.2017.06.082](#)
13. Breast Cancer Res., 19, 20, (2017), [DOI: 10.1186/s13058-017-0813-x](#)
14. Appl Microbiol Biotechnol., 101, 7371, (2017), [DOI: 10.1007/s00253-017-8471-6](#)
15. J Control Release., 265, 66, (2017), [DOI: 10.1016/j.jconrel.2017.04.027](#)
16. J Mol Biol. (1975), [DOI: 10.1016/0022-2836\(75\)90231-4](#)

### Product attributes

CAS number	58880-05-0
Excitation/Emission	510/600 nm (after photocrosslinking to nucleic acid)
Apoptosis/viability marker	Dead cell stain
Assay type/options	Live/dead discrimination, Viability PCR
Detection method/readout	PCR/qPCR
Cell permeability	Membrane impermeant
Molecular weight	420
Storage Conditions	Store at 2 to 8 °C or below, Protect from light, Desiccate