



Glowing products for science

Call us : [800-304-5357](tel:800-304-5357)

## Mix-n-Stain™ FITC Antibody Labeling Kit



Label your antibody with fluorescein (FITC) just 15 minutes without a purification step. Labeling tolerates many common buffer components including BSA and ascites.

### Product attributes

Excitation/Emission	494/518 nm
---------------------	------------

### Product Description

Mix-n-Stain™ FITC Antibody Labeling Kits allow you to label  $\leq 5$  ug up to 100 ug of your antibody with fluorescein (FITC) in just 15 minutes, with minimal hands-on time and no purification.

#### Unrivaled Convenience

Mix-n-Stain™ labeling tolerates common antibody buffer formulations. An ultrafiltration centrifuge vial is provided to quickly remove interfering substances like glycerol, or to concentrate your antibody. With a slight modification in protocol, antibodies can be labeled in the presence of BSA, gelatin, or ascites fluid. The presence of other proteins like BSA or gelatin in the labeling reaction has minimal effect on background fluorescence, because any labeled non-antibody proteins readily wash away during immunofluorescence staining.

#### Choose the Right Labeling Kit for Your Antibody

This kit is for labeling an antibody with the traditional green fluorescent dye fluorescein (FITC). For a green dye with superior brightness, photostability, and insensitivity to pH, see our [CF@488A Mix-n-Stain™ Antibody Labeling Kit](#). Learn more about our next-generation [CF@ Dyes](#).

Mix-n-Stain™ Antibody Labeling Kits are very simple to use (see the workflow overview above). But before you begin, you must check that your antibody meets the compatibility requirements for labeling, and choose the right labeling protocol. Download the updated and expanded [Product Protocol](#) for a labeling check list to help you select the right kit size and labeling protocol for your antibody. Also see our [Kit Compatibility and Protocol Selection Flowchart](#).

We offer Mix-n-Stain™ Kits for labeling antibodies, nanobodies, and small ligands with a wide selection of fluorescent dyes, enzymes, and other labels. See our full selection of Mix-n-Stain™ Antibody Labeling Kits below.

# Antibody & Protein Labeling Kits

Product Name	Label Options	Labeling Scale	Labeling Time	Features
<a href="#">Mix-n-Stain™ CF® Dye Antibody Labeling Kits</a>	CF® Dyes	≤5-20 ug IgG 20-50 ug IgG 50-100 ug IgG	~ 15 min.	<ul style="list-style-type: none"> <li>• Rapid, simple labeling</li> <li>• No purification</li> <li>• Tolerates BSA &amp; other additives</li> <li>• Dyes for super-resolution, spectral flow, &amp; NIR detection</li> </ul>
<a href="#">Mix-n-Stain™ FITC Antibody Labeling Kits</a>	FITC			
<a href="#">Mix-n-Stain™ Cyanine Dye Antibody Labeling Kits</a>	Cyanine 555 Cyanine 647			
<a href="#">Mix-n-Stain™ Biotin Antibody Labeling Kits</a>	Biotin			
<a href="#">Mix-n-Stain™ Digoxigenin Antibody Labeling Kits</a>	Digoxigenin (DIG)			
<a href="#">Mix-n-Stain™ DNP Antibody Labeling Kits</a>	DNP			
<a href="#">Mix-n-Stain™ Maxi Antibody Labeling Kits</a>	CF® Dyes Cyanine Dyes	1 mg IgG	~ 30 min.	
<a href="#">Mix-n-Stain™ CF® Dye IgM Antibody Labeling Kits</a>	CF® Dyes FITC	25 ug IgM 100 ug IgM	15-30 min.	<ul style="list-style-type: none"> <li>• Rapid, simple labeling</li> <li>• No purification</li> </ul>
<a href="#">Mix-n-Stain™ STORM CF® Dye Antibody Labeling Kits</a>	CF® Dyes	50 ug IgG	~ 30 min.	<ul style="list-style-type: none"> <li>• Optimized to provide low 1-2.5 DOL</li> <li>• No purification</li> <li>• Optimized for Nanobodies®</li> </ul>
<a href="#">Mix-n-Stain™ Nanobody Labeling Kits</a>	CF® Dyes	5-20 ug Nanobody® 20-50 ug Nanobody®	~ 30 min.	<ul style="list-style-type: none"> <li>• Tolerates BSA &amp; other additives</li> <li>• Label SNAP®, CLIP™, HaloTag® &amp; other small ligands</li> </ul>
<a href="#">Mix-n-Stain™ Small Ligand Labeling Kits</a>	CF® Dyes	0.1 umol small ligand	~ 30 min.	<ul style="list-style-type: none"> <li>• Dye options for surface or intracellular targets</li> </ul>
<a href="#">VivoBrite™ Rapid Antibody Labeling Kits for Small Animal In Vivo Imaging</a>	Near-IR CF® Dyes	3 x 1 mg IgG	~ 2 hours	<ul style="list-style-type: none"> <li>• Dyes, buffers, &amp; spin vials for labeling + purification</li> <li>• Sterile syringes, filters and storage vials for in vivo use</li> </ul>
<a href="#">CF® Dye &amp; Biotin SE Protein Labeling Kits</a>	CF® Dyes Biotin	3 x 1 mg protein	~ 2 hours	<ul style="list-style-type: none"> <li>• Dyes, buffers, &amp; spin vials for labeling + purification</li> </ul>

## References

Download a list of [Mix-n-Stain™ references](#).

This datasheet was generated on May 9, 2026 at 02:15:09 AM. Visit product page to check for updated information before use.

Product link: <https://biotium-woo.supremeclients.com/product/mix-n-stain-fitc-antibody-labeling-kit/>