



Glowing products for science

Call us : [800-304-5357](tel:800-304-5357)

BactoView™ Dead Stains

Highly-selective dead cell stains for bacteria cultures. The stains are suitable for both gram-positive and gram-negative strains.



Product Description

BactoView™ Dead Stains are novel DNA binding dyes for live/dead discrimination in bacteria. These bright fluorogenic DNA binding dyes are cell membrane-impermeant, so they selectively stain dead bacteria with compromised cell membranes without the need to wash after staining.

- Highly-selective dead cell stains for bacteria cultures
- Superior live/dead differentiation for gram-positive strains than traditional dyes
- Also provide excellent results in gram-negative strains
- Quick 30-minute staining, no wash required
- Available in 7 colors from green to near-IR

Note: As of 6/16/2025, the formulation of BactoView™ Dead Stains (Cat. No. 40107-40113) was updated from 500X in water to 1000X in DMSO.

Highly-Selective Dead Cell Staining for Gram-Positive and Gram-Negative Strains

Traditional vital nucleic acid dyes like propidium iodide or ethidium homodimer are efficiently excluded from live gram-negative bacteria but are taken up by live gram-positive bacteria, resulting in high background in live cells and poor live/dead discrimination. BactoView™ Dead Stains have novel chemical structures that are efficiently excluded from both gram-positive and gram-negative strains, for highly selective live/dead discrimination. In addition to staining dead bacteria, BactoView™ Dead Stains also stain *Bacillus subtilis* endospores, but with dimmer fluorescence. BactoView™ Dead Stains have low fluorescence until they bind DNA, allowing bright, no-wash staining. The stains are available with a wide selection of emission wavelengths ranging from green to near-infrared, for microscopy or flow cytometry analysis.

Product attributes

| | |
|-----------------------------------|--|
| Apoptosis/viability marker | Dead cell stain |
| For live or fixed cells | For live/intact cells |
| Detection method/readout | Fluorescence microscopy, Flow cytometry |
| Assay type/options | Endpoint assay |
| Colors | Green, Red, Far-red, Near-infrared |
| Storage Conditions | Store at -10 to -35 °C, Protect from light |

BactoView™ Dead Stains & Viability Kits

| Product Name | Ex/Em (nm) | Detection Channel | Size (1000X in DMSO) | Catalog No. |
|--|-----------------------------------|--------------------------------------|----------------------|-------------------------|
| BactoView™ Dead 500/515 100 uL | 497/515 40107 | FITC | 20 uL | 40107-T |
| BactoView™ Dead 560/570 100 uL | 559/570 40108 | Rhodamine, PI, PE | 20 uL | 40108-T |
| BactoView™ Dead 570/585 100 uL | 572/583 40109 | Rhodamine, PI, PE | 20 uL | 40109-T |
| BactoView™ Dead 600/615 100 uL | 603/613 40110 | Texas Red® or PE-Texas Red® | 20 uL | 40110-T |
| BactoView™ Dead 655/670 100 uL | 653/671 40111 | Cy®5, APC | 20 uL | 40111-T |
| BactoView™ Dead 690/710 100 uL | 683/707 40112 | Cy®5.5 | 20 uL | 40112-T |
| BactoView™ Dead 760/780 100 uL | 759/780 40113 | Cy®7, APC-Cy®7 | 20 uL | 40113-T |
| BactoView™ Viability Kit (Green/Red) | Green (498/522)/Red (572/583) | FITC (Green)/Rhodamine, PI, PE (Red) | 1 kit | 32019 |
| BactoView™ Viability Kit (Green/Far-Red) | Green (498/522)/Far-Red (653/671) | FITC (Green)/Cy®5, APC (Far-Red) | 1 kit | 32020 |

Texas Red is a registered trademark of Thermo Fisher Scientific; Cy Dye is a registered trademark of Cytiva.

Staining of Gram-Negative *E. Coli*

Staining of Gram-Positive *B. Subtilis*

BactoView™ Dead Also Stains *B. Subtilis* Endospores

Biotium also offers [BactoView™ Viability Kits](#), which include a choice of red or far-red BactoView™ Dead Stain for dead bacteria and BactoView™ Viability Green Counterstain to stain all bacteria. BactoView™ Dead Stains also can be combined with fluorescent Gram stains like our [CF® Dye WGA Conjugates](#).

We also offer [BactoView™ Live Green](#) and [BactoView™ Live Red](#) for staining live bacteria. Note that BactoView™ Stains cannot be used to distinguish bacteria from eukaryotic cells, because they will stain other cell types as well. For bright and optimized labeling of bacterial endospores, see our [BactoSpore™ Bacterial Stains](#).

For staining mammalian cells, see our [NucSpot® Nuclear Stains](#) for live/dead discrimination or nuclear counterstaining of fixed mammalian cells. For live nuclear staining of mammalian cells, see our [NucSpot® Live Stains](#). Also, view our [Cellular Stains Table](#) for more information on how our dyes stain various organisms.

This datasheet was generated on May 8, 2026 at 01:15:16 AM. Visit product page to check for updated information before use.

Product link: <https://biotium-woo.supremeclients.com/product/bactoview-dead-stains/>