

## TPEN

TPEN (Tetrakis-(2-pyridylmethyl)ethylenediamine) is an excellent heavy metal chelator while it does not affect calcium, magnesium, sodium, or potassium concentrations. Thus, TPEN is a useful tool to identify effects of heavy metals (zinc, iron, copper, manganese, etc.) on fluorescent indicators for other ions.



### Product attributes

CAS number	16858-02-9
Cell permeability	Membrane permeant

## Product Description

TPEN (Tetrakis-(2-pyridylmethyl)ethylenediamine) is an excellent heavy metal chelator while it does not affect  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$ ,  $\text{Na}^{+}$ , or  $\text{K}^{+}$  concentrations. Thus, TPEN is a useful tool to identify effects of heavy metals ( $\text{Zn}^{2+}$ ,  $\text{Fe}^{2+}/\text{Fe}^{3+}$ ,  $\text{Cu}^{2+}$ , and  $\text{Mn}^{2+}$ , etc.) on fluorescent indicators of  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$ ,  $\text{Na}^{+}$ , or  $\text{K}^{+}$ .

- Off-white solid soluble in ethanol, DMSO or DMF
- Store at 4 °C
- $\text{C}_{26}\text{H}_{28}\text{N}_6$
- MW: 425
- [16858-02-9]

As the indicator does not covalently bind to cellular components, it may be actively effluxed from the cell by organic anion transporters. The rate of efflux increases with temperature, and may vary between cell types, resulting in variable retention times of a few minutes to hours. Experiments using indicators in cells usually are performed within one or two hours of loading, but it may be possible to re-load cells with indicator if needed. The organic anion transporter inhibitor [Probenecid \(#50027\)](#) can be used to slow the rate of indicator efflux from cells.

## References

1. J Biol Chem, 260, 2719 (1985), [PMID: 3919006](#)
2. Methods Cell Biol, 99, 113, (2021), [DOI: 10.1016/B978-0-12-374841-6.00005-0](#)

This datasheet was generated on May 8, 2026 at 06:47:54 AM. Visit product page to check for updated information before use.  
Product link: <https://biotium-woo.supremeclients.com/product/tpen-full-name-tetrakis-2-pyridylmethylethylenediamine/>