

Biotin-X cadaverine, trifluoroacetate salt

Biotin-X cadaverine trifluoroacetic acid salt (5-(((N-(biotinoyl)amino)hexanoyl)amino) pentylamine, trifluoroacetic acid salt) can be coupled to activated carboxylic acid esters or sulfonyl chlorides in the presence of a base or in a basic buffer. Also see biotin-X cadaverine free base ([90080](#)).




Product attributes

Product Description

Biotin-X cadaverine trifluoroacetic acid salt (5-(((N-(biotinoyl)amino)hexanoyl)amino) pentylamine, trifluoroacetic acid salt) can be coupled to activated carboxylic acid esters or sulfonyl chlorides in the presence of a base or in a basic buffer. Also see biotin-X cadaverine free base ([90080](#)).

- Tan solid soluble in DMSO
- Store at -20°C or 4°C
- C₂₃H₄₀F₃N₅O₅S
- MW: 555.65



References

1. J Histochem Cytochem 38, 377 (1990).
2. Ann NY Acad Sci 463, 214 (1984).
3. Biochem J 251, 935 (1988).

1. J Histochem Cytochem 38, 377 (1990).
2. Ann NY Acad Sci 463, 214 (1984).
3. Biochem J 251,935 (1988).

This datasheet was generated on May 8, 2026 at 04:31:42 AM. Visit product page to check for updated information before use.

Product link: <https://biotium-woo.supremeclients.com/product/biotin-x-cadaverine-trifluoroacetate-salt-5-n-biotinoylaminohexanoylamino-pentylamine-trifluoroacetic-acid-salt/>