

## Mag-Fluo-4, AM Ester

Membrane-permeant AM ester form of Mag-Fluo-4, a fluorescent magnesium and calcium ion indicator that is an analog of Fluo-4.



### Product attributes


Cell permeability	Membrane permeant
Indicator type	Non-ratiometric
Excitation/Emission	493/517 nm (high Mg <sup>2+</sup> or Ca <sup>2+</sup> )

## Product Description

Membrane-permeant AM ester form of Mag-Fluo-4, a fluorescent ion indicator and analog of Fluo-4.

Because of the relatively low water solubility of the AM ester, the mild detergent [Pluronic® F-127 \(cat# 59004\)](#) is often used as a dispersing agent to facilitate loading. The indicator has an affinity for magnesium (K<sub>d</sub> = 4.7 mM) and low-affinity for calcium (K<sub>d</sub> = 22 uM) with absorbance/emission at 493/517 nm (high Mg<sup>2+</sup>).

- λ<sub>Ex</sub>/λ<sub>Em</sub> (high Mg<sup>2+</sup> or Ca<sup>2+</sup>) (after hydrolysis): 493/517 nm
  - Orange solid soluble in DMSO
  - Store -20°C, desiccate and protect from light
  - C<sub>37</sub>H<sub>33</sub>F<sub>2</sub>NO<sub>18</sub>
  - MW: 818



BAPTA-based ion indicators like Mag-Fluo-4 have been shown to be fixable in situ by [EDC/EDAC \(cat# 59002\)](#). The fixation of indicator dyes is useful for downstream immunofluorescence and IHC studies ([Cell Calcium 1997, 21\(3\), 175](#)).

As the indicator does not covalently bind to cellular components, it may be actively effluxed from the cell by organic anion transporters. The rate of efflux increases with temperature, and may vary between cell types, resulting in variable retention times of a few minutes to hours. Experiments using indicators in cells usually are performed within one or two hours of loading, but it may be possible to re-load cells with indicator if needed. The organic anion transporter inhibitor [Probenecid \(#50027\)](#) can be used to slow the rate of indicator efflux from cells.

Pluronic is a registered trademark of BASF.

## References

1. Biophys J. (2009) Oct 7; 97(7):1864 [doi: 10.1016/j.bpj.2009.07.021](https://doi.org/10.1016/j.bpj.2009.07.021)
2. Biochem J. (2002) Oct 1; 367(Pt1):137 [doi: 10.1042/BJ20020305](https://doi.org/10.1042/BJ20020305)
3. Biochimica et Biophys Acta - Mol Cell Res. (2017) June; 1864(6):977 [doi.org/10.1016/j.bbamcr.2016.11.026](https://doi.org/10.1016/j.bbamcr.2016.11.026)
4. Biochem & Biophys Res Comm. (2014) November 28; 454(4):572 [doi.org/10.1016/j.bbrc.2014.10.125](https://doi.org/10.1016/j.bbrc.2014.10.125)
5. Methods Cell Biol, 99, 113, (2021), [DOI: 10.1016/B978-0-12-374841-6.00005-0](https://doi.org/10.1016/B978-0-12-374841-6.00005-0)