

## GFAP Recombinant Alpaca VHH (SdAb2409.GFAP) - MiniMab™



VHH nanobody targeting glial fibrillary acidic protein (GFAP), part of our MiniMab™ SdAb series that have been engineered for optimal conjugate performance.

### Product Description

GFAP Recombinant Alpaca VHH (SdAb2409.GFAP) recognizes the glial fibrillary acidic protein (GFAP) expressed in neural tissues. This high-affinity single-domain antibody (SdAb), also known as camelid VHH or Nanobody®, is part of our MiniMab™ series of highly optimized conjugated probes. The nanobody has been validated for immunofluorescence microscopy and is available conjugated to CF® Dyes.

#### Features of MiniMab™ Nanobodies single-domain antibodies

- Superior to conventional antibodies: deeper tissue penetration, higher solubility and stability, and faster staining
- Minimal epitope-dye displacement—perfect for super-resolution imaging
- Specifically developed and optimized for immunofluorescence
- Labeled with bright, photostable CF® Dyes, including near-infrared CF®740 option
- Available as conjugates with Biotium's best-in-class dyes for STORM

GFAP is a 49-kDa type III intermediate filament protein found in neural tissues, serving as a marker that differentiates astrocytes from other glial cells during the development of the central nervous system. Although three splice variants of GFAP have been identified, the alpha isoform is the most abundantly expressed in astrocytes. GFAP is capable of co-assembling with vimentin and nestin in astrocytes, though these interactions are not essential for its filament formation. Similar to other intermediate filaments, GFAP assembly is regulated by phosphorylation and dephosphorylation of its N-terminal domain. Mutations in the GFAP gene have been linked to Alexander disease, and elevated GFAP expression has been observed in certain tumors derived from glial cells.

Learn more about [CF® Dyes for super-resolution](#) as well as our innovative reagents for [immunofluorescence microscopy](#); this includes our [NucSpot® Nuclear Stains](#) for bright and nuclear-specific staining in a wide color selection, and [Cytoliner™ Fixed Cell Membrane Stains](#) for robust membrane staining in formaldehyde-fixed cells.

[View our full selection of primary and secondary antibodies](#) available with bright CF® Dyes and other labels.

### Product attributes

Antibody number	N001
Antibody type	MiniMab™ SdAb (VHH)
Clonality	Recombinant single-domain antibody
Host species	Alpaca
Clone	SdAb2409.GFAP
Antibody reactivity (target)	GFAP
Synonyms	Astrocyte or Intermediate Filament Protein, Glial Fibrillary Acidic Protein (GFAP)
Species reactivity	Human, Mouse, Rat
Human gene symbol	GFAP
Entrez gene ID	14580
Molecular weight	50kDa
Cell/tissue expression	Glia
Verified antibody applications	IF (verified)
Positive control	Brain, Retina
Antibody application notes	Immunofluorescence (HeLa): 1-5 ug/mL Optimal concentration to be determined by end-user.
Antibody research areas	Neuroscience
Antibody/conjugate formulation	Conjugates: PBS/0.1% rBSA/0.05% azide
Shelf life	Guaranteed for at least 24 months from date of receipt when stored as recommended
Storage Conditions	Store conjugates at 2 °C to 8 °C, Protect fluorescent conjugates from light
Shipping condition	Room temperature
Regulatory status	For research use only (RUO)
Product origin	Recombinant alpaca VHH produced in E.coli, Recombinant BSA produced in Chinese hamster ovary cells

## GFAP Recombinant Alpaca VHH (SdAb2409.GFAP) – MiniMab™

Conjugation	Ex/Em	Conc.	STORM <sup>1</sup> compatibility	Catalog No.	Dye Features
<a href="#">CF®488A</a>	490/516 nm	100 ug/mL	Yes	<a href="#">N001-488A-200UL</a>	<a href="#">CF®488A Features</a>
<a href="#">CF®498</a>	498/519 nm	100 ug/mL	Yes	<a href="#">N001-498-200UL</a>	
<a href="#">CF®568</a>	562/584 nm	100 ug/mL	Yes	<a href="#">N001-568-200UL</a>	<a href="#">CF®568 Features</a>
<a href="#">CF®583R</a>	585/609 nm	100 ug/mL	Yes	<a href="#">N001-583R-200UL</a>	<a href="#">CF®583R Features</a>
<a href="#">CF®647</a>	652/668 nm	100 ug/mL	Yes	<a href="#">N001-647-200UL</a>	<a href="#">CF®647 Features</a>
<a href="#">CF®660C</a>	667/685 nm	100 ug/mL	Yes	<a href="#">N001-660C-200UL</a>	<a href="#">CF®660C Features</a>
<a href="#">CF®680</a>	681/698 nm	100 ug/mL	Yes	<a href="#">N001-680-200UL</a>	<a href="#">CF®680 Features</a>
<a href="#">CF®740</a>	742/767 nm	100 ug/mL	No	<a href="#">N001-740-200UL</a>	<a href="#">CF®740 Features</a>

<sup>1</sup> STORM: Stochastic optical reconstruction microscopy. [Learn more about CF® Dyes for super-resolution.](#)  
 NANOBODY is a registered trademark of ABLYNX.