

## Bovine Serum Albumin Fraction V

Can be used as a stabilizing agent for proteins and enzymes, or added to blocking buffers to minimize non-specific antibody binding in immuno-detection procedures such as western blotting and cell staining.



### Product attributes

Storage Conditions
Store at 2 to 8 °C

## Product Description

Bovine serum albumin (BSA) Fraction V has uses in many applications such as stabilizing proteins in solution, blocking non-specific binding in immuno-assays, and as a protein standard in various protein assays.

- Stabilizer for proteins and enzymes in solution
- Blocking agent for immuno-assays (ie, western blot, IHC, IF)
- Minimizes non-specific antibody binding
- IgG- and protease-free
- 50 g of solid powder

Bovine serum albumin (BSA) is commonly used as a stabilizing agent for proteins and enzymes, including dilute solutions of antibody. It is also used as a blocking agent to reduce non-specific antibody binding in immuno-detection procedures such as western blotting, immunofluorescence, and IHC. It is also suitable for use in protease sensitive assays such as RIA and EIA, as well as nucleic acid hybridization. This BSA is IgG- and protease-free.

This item is only available for customers in the USA.

## References

ACS Omega (2020) Feb 26;5(9):4673-4681 [doi: 10.1021/acsomega.0c00115](https://doi.org/10.1021/acsomega.0c00115).

This datasheet was generated on May 9, 2026 at 11:04:38 AM. Visit product page to check for updated information before use.  
Product link: <https://biotium-woo.supremeclients.com/product/bovine-albumin-serum-fraction-v/>