

Biotin-11-dUTP, 1 mM Solution

Biotin-11-dUTP can be enzymatically incorporated into DNA via nick translation, random priming, or 3'-end terminal labeling. The number '11' is the number of atoms in the linker between biotin and dUTP.




Product attributes

CAS number 86303-25-5

Product Description

Biotin-11-dUTP can be enzymatically incorporated into DNA via nick translation, random priming, or 3'-end terminal labeling. The number '11' is the number of atoms in the linker between biotin and dUTP. Biotium also offers biotin-16-dUTP ([40022](#)) and biotin-20-dUTP ([40030](#)). The length of the linker affects the incorporation efficiency of the biotin-dUTP probe into DNA using DNA polymerases, and it also affects biotin/avidin or biotin/streptavidin. In general, the shorter the linker, the more efficiently the biotin-dUTP is incorporated into DNA by DNA polymerases. On the other hand, the longer the linker, the better biotin can interact with avidin or streptavidin. Biotin-11-dUTP is also available as a lyophilized solid ([40029-1](#)).

- 1 mM solution in pH 7.5 Tris-HCl buffer
- Store at -20 °C
- C₂₈H₄₁N₆O₁₇P₃S Li₄
- MW: 886.5
- [86303-25-5]



This datasheet was generated on May 8, 2026 at 07:49:17 PM. Visit product page to check for updated information before use.
 Product link: <https://biotium-woo.supremeclients.com/product/biotin-11-dutp-1-mm-in-ph-7-5-tris-hcl-buffer/>