

## Aminoxy-biotin

Aminoxy-biotin can be used to covalently attach biotin to aldehyde or ketone groups on polysaccharides, glycoproteins or antibodies. Aminoxy-biotin has also been used to efficiently label cell-surface glycans on living animal cells.




### Product attributes

Chemical reactivity (reacts with)	Aldehydes/ketones
Functional group	Aminoxy (hydroxylamine)
Storage Conditions	Store at 2 to 8 °C

## Product Description

Aminoxy (aka hydroxylamine or aminoxyacetamide) biotin readily reacts with aldehyde or ketone groups to form a stable oxime linkage under mild conditions, without the use of reducing agents.

- Stable labeling of aldehyde or ketone groups on polysaccharides, glycoproteins or antibodies.
- Mild labeling conditions, at neutral to acidic pH, in aqueous solutions.
- Label cell-surface sialic acid groups on glycans on living animal cells.
- White to off-white solid soluble in DMSO.
- MW: 515.55.



We also offer [CF@ Dye Aminoxy](#) derivatives that are excellent options for fluorescent labeling. [CF@ dyes](#) are Biotium's line of next-generation fluorescent dyes that have improved brightness, photostability and water solubility compared to other fluorescent dyes.

### Advantages of Aminoxy Labeling

Aminoxy groups react with molecules containing aldehyde or ketone groups to form a stable oxime linkage in aqueous solutions, at neutral/acidic pH. The reaction is fast and can be further accelerated using a catalyst like aniline. Aminoxy reagents are therefore superior to hydrazides, which also react with aldehydes or ketones but form unstable hydrazone linkages. Aminoxy reagents offer a convenient and rapid way to label glycoproteins with a detectable tag, such as a fluorescent dye or biotin, under mild conditions. The labeling process involves two steps, first introducing aldehyde groups into the glycoproteins by mild periodate oxidation, followed by treating the functionalized proteins with an aminoxy reagent. Aminoxy labeling of antibody glycosylation sites can be used as an alternative to succinimidyl ester labeling of amines, for antibodies where amine labeling affects the antibody binding affinity. Aminoxy biotin has also been used to efficiently label cell-surface sialic acid-containing glycans on living animal cells.

## References

1. Sci Rep (2019) 9, 9018. [DOI:10.1038/s41598-019-44957-y](https://doi.org/10.1038/s41598-019-44957-y)
2. PNAS (2000) 97, 686-691. [DOI: 10.1073/pnas.97.2.686](https://doi.org/10.1073/pnas.97.2.686)
3. Biochemistry (1993) 32, 8276-8283. [DOI:10.1021/bi00083a031](https://doi.org/10.1021/bi00083a031)
4. Biochemistry (1992) 31, 3703-3708. [DOI: 10.1021/bi00129a020](https://doi.org/10.1021/bi00129a020)

This datasheet was generated on May 8, 2026 at 07:19:27 AM. Visit product page to check for updated information before use.  
Product link: <https://biotium-woo.supremeclients.com/product/aminoxy-biotin/>