Glowing products for science

Call us : 800-304-5357

Flow Cytometry Fixation/Permeabilization Kit

A kit containing fixation and permeabilization buffers used in staining intracellular antigens for flow cytometry.



Product attributes

Storage Conditions	Store at room temperature

Product Description

The Flow Cytometry Fixation/Permeabilization Kit contains optimally formulated buffers for fixation and permeabilization of cells for immunofluorescence staining of intracellular antigens for analysis by flow cytometry.

Note: the fixation buffer in this kit is prepared from methanol-stabilized formaldehyde. For methanol-free fixative, see our Paraformaldehyde, 4% in PBS, Ready-to-Use Fixative (22023).

See our Flow Cytometry Technology Page to see all of our products for flow cytometry.

This datasheet was generated on November 8, 2026 at 10:50:53 AM. Visit product page to check for updated information before use.
Product link: https://biotium-woo.supremeclients.com/product/flow-cytometry-fixationpermeabilization-kit/